

Alternating Low Air Loss Mattress Replacement



Matrix ALAL

- Twenty 8" High Nylon Polyurethane Low Air Loss Cells, Featuring an AVB
- Alternating System, with 4 Static head Zone Cells
- Fully inflates mattress in under two minutes
- Programmable 5/10/15/20 minute Alternation Cycles at the Push of a Button for Easy Use
- Patient comfort control panel features an adjustable Soft/Firm key to provide variable pressure to accommodate the patient's weight or comfort up to 500 pounds
- Breathable Soft, 4-Way Stretch Top Cover, Fully Removable
- 2" Vinyl Foam Base Provides for Extra Comfort and additional protection while patient is transported through facility
- Custom Internal Manifold System with "Quick Disconnect" Feature for Optimum Safety
- Lockout Feature
- Surface has 8" minimum depth
- Bariatric sizes available up to 48
- Bariatric weight capacity up to 1,000 pounds

CUSTOM MEDICAL
SOLUTIONS

(8 6 6) 3 5 0 - 5 6 4 0

www.custommedicalsolutions.com